

Instructions for the shipment of samples

1. SUGGESTED TYPE OF CONTAINERS FOR SAMPLES

Sender must send the samples in sealed containers, unbreakable and suitable for various types of wastes, either solid, liquid, mud or debris. It is suggested to use the following containers:

• Liquids: 1 liter / 30 oz polyethylene bottle

Solids or semi-solids:
 Debris:
 500 ml / 15 oz polyethylene wide-mouth jar
 20" x 3" hi-density nylon bag with tie

Samples must be placed in a cardboard box for transportation.

2. LABELLING

The box must be labelled with the following information:

- Name and address of generator
- Description of waste
- Number of samples and individual weight
- Date
- Destination : Stablex Canada, c/o Pierre-Olivier Gagné

SHIPMENT

From	Sample Gross weight		Sent to:	
Canada	Equal or less than 10 kg or 10 liters	Stablex Canada Inc c/o Pierre-Olivier Gagné 760 Industrial Blvd. Blainville QC J7C 3V4		
USA	Equal or less than 5 kg or 5 liters	Stablex Canada Inc c/o Pierre-Olivier Gagné 760 Industrial Blvd. Blainville QC J7C 3V4 (450) 430-9230	OR	EQ Northeast * c/o Richard Bellefeuille 10 Terminal Drive Latham NY 12110 (508) 243-4846
Via FedEx or UPS ground services				

Ref: http://www.fedex.com/downloads/hk_english/packagingtips/pointers.pdf

4. QUESTIONS

Any question regarding instructions for the shipment of samples can be directed to Pierre-Olivier Gagné, at pierreolivier.gagne@stablex.com or 1-800 782-2539 ext. 884744.

SCI-1000-45/dec 2022 Reference: SCI-20-205

^{*} Sender must advise Stablex Customer Service before shipping samples to EQ Northeast.